# Data Quality Monitoring for High Energy Physics (DQM4HEP) Version 03-02-00

# R. Été, A. Pingault, L. Mirabito

Université Claude Bernard Lyon 1 - Institut de Physique Nucléaire de Lyon / Ghent University

26 avril 2016







Architecture and API

Conclusion and plans

DQM4HEP: an online monitoring system for data quality

### Key points

- Event distributed system : server/client paradigm
- Set of interfaces for data analysis, adapted to DQM purpose
- Histogram distributed system
- Visualization interface (Qt GUI)
- Large scale remote process management
- Generic IO support for any edm (opt. LCIO)
- Full size HEP experiment to single detector prototype design
- ELog interface

Set of interfaces inspired from CMS DQM system (monitor elements, collectors).

Application flow inspired from ALICE DQM system, AMORE (cycles).



DQM4HEP packages

One location: https://github.com/DQM4HEP

Webpage: dqm4hep.github.io

# The main package: DQM4HEP

Installation package for sub-packages (CMake).

Sub-packages:

- dim: Distributed Information Management (Delphi). Manage client/server communications
- dimjc: DIM Job Control (L. Mirabito). Remote process management using dim.
- jsoncpp : Json I/O for dimjc
- streamlog : logging library (used in ILCSOFT)
- DQMCore: Core part of the DQM system. Client/server interfaces, analysis, IO, run control interface, plugin management ...
- DQMViz: Qt visualization interfaces. Job control gui client, monitoring gui client, run control server gui (standalone).
- LCIO: Linear Collider IO. Build support for LCIO streamer

### Forseen packages:

- xdrstream : Generic Xdr serializer
- xdrlcio: Lcio serialization using xdrstream (buffer -> socket)
- DQM4ILC : ILC specific implementation (detector prototypes modules, marlin helper, ...)

Installation

#### Installation mode

Designed to be built standalone or using ILCSOFT.

Basic install requires ROOT.

Full install with DQMViz requires Qt and ROOT compiled with -enable-qt option.


### Standalone mode:

- Basic install : dim, dimjc, jsoncpp, streamlog, DQMCore
- Full install : + DQMViz, LCIO


### ILCSOFT mode:

- Basic install : dim, dimjc, jsoncpp, DQMCore
- Full install : + DQMViz

Global workflow



DAQ interface example : SDHCAL DAQ interface



Module applications - analysis module

#### Purpose

- Receive events from a collector server and process them
- Produce monitor elements (histograms, scalars, generic TObject)
- Follow the run control signals (SOR, EOR)
- Init : Initialize the application : load dlls, declare services, etc ... Wait for a SOR
- Start of run: start cycles loop, open archive
- Start of cycle : start a cycle of 'process event'
- Process event : Process incoming event, fill monitor elements, etc ...
- End of cycle: send subscribed monitor elements, update archive (opt).
- End of run : Wait for SOR, close archive (opt).
- End: Clean and exit module.


To implement online DQM analysis, user must implement the DQMAnalysisModule interface. A shared library must be build and loaded in the application using the plugin system (export DQM4HEP\_PLUGIN\_DLL=libMyModule.so).

Use dqm4hep\_start\_analysis\_module to start an analysis module.



Analysis module application flow

Module API



Module API

### Monitor element management

- Memory management
- Booking methods (book, delete, from xml)
- Access management via an internal map ("name" -> element)
- Directory management : mkdir, cd, ls, rmdir, pwd

### Quality tests

Quality tests used to **evaluate the quality** of a particular monitor element.

Example of quality tests:

- Kolmogorov test with a reference histogram
- χ<sup>2</sup> test after a histogram fit

Quality tests can also be user defined (via DQMQualityTestFactory class).

Quality tests results (DQMQualityTestResult class) stored within monitor element and both sent to the collector.

The module API provides functions for :

- User QTests registration
- QTests assignment to monitor elements (also from xml)
- QTest results access (read only)

Analysis module - Example

```
// ExampleModule.h
class ExampleModule : public DQMAnalysisModule
public:
  ExampleModule(); // requiered by plugin system
 ~ExampleModule():
  StatusCode initModule();
  StatusCode readSettings(const TiXmlHandle &
        xmlHandle):
  StatusCode startOfRun(DQMRun *const pRun);
  StatusCode startOfCvcle():
  StatusCode processEvent(DQMEvent *const pEvent):
  StatusCode endOfCycle();
  StatusCode endOfRun(DQMRun *const pRun);
  StatusCode endModule():
private
  DOMMonitorFlement
                      *m pNHitElement:
```

```
// ExampleModule.cc
#include "ExampleModule.h"
// declare plugin in the system
DQM PLUGIN DECL( ExampleModule, "ExampleModule")
// create and enter dir. Book histogram
StatusCode ExampleModule::initModule()
  DQMModuleApi::mkdir(this, "/Histograms");
  DQMModuleApi::cd(this. "/Histograms"):
  DQMModuleApi::bookIntHistogram1D(this,
        m pNHitElement, "NHit",
    "Number of hits", 1200, 0, 1199);
  return STATUS CODE SUCCESS;
// get Icio event and fill your histogram !
StatusCode ExampleModule::processEvent(DQMEvent *
      const pEvent)
  EVENT::LCEvent *pLCEvent =
    pEvent->getEvent<EVENT::LCEvent>();
  if (!pLCEvent)
    return STATUS CODE FAILURE:
    // get number of hits from a collection
  m pNHitElement->get<TH1I>()->Fill(NHit);
  return STATUS CODE SUCCESS:
```

Gui visualisation

# Gui interfaces for DQM client developed :

- Run control, job control, online monitoring
- Written with Qt4 framework
- Easily configurable with json and xml.

#### Run Control GUI



 Parametrisation of run with run number, detector name, run description and parameters

#### Run Control GUI



- Parametrisation of run with run number, detector name, run description and parameters
- Send SOR and EOR signals

#### Run Control GUI




- Parametrisation of run with run number, detector name, run description and parameters
- Send SOR and EOR signals
- Control run status (State, Started/Stopped time)

#### Run Control GUI




- Parametrisation of run with run number, detector name, run description and parameters
- Send SOR and EOR signals
- Control run status (State, Started/Stopped time)
- Every action is logged for easy information overview

Job Control GUI




Load and display a list of applications (Collectors, Modules, etc.) available on different hosts

Job Control GUI




- Load and display a list of applications (Collectors, Modules, etc.) available on different hosts
- Displays informations(Name, Host, PID, Status, etc.) about applications

Job Control GUI



- Load and display a list of applications (Collectors, Modules, etc.) available on different hosts
- Displays informations(Name, Host, PID, Status, etc.) about applications
- Infos can be updated in "real time"

Job Control GUI




- Load and display a list of applications (Collectors, Modules, etc.) available on different hosts
- Displays informations(Name, Host, PID, Status, etc.) about applications
- Infos can be updated in "real time"
- Manage Applications (Start/Kill/Restart) with contextual menu



14 / 18

Job Control GUI




- Load and display a list of applications (Collectors, Modules, etc.) available on different hosts
- Displays informations(Name, Host, PID, Status, etc.) about applications
- Infos can be updated in "real time"
- Manage Applications (Start/Kill/Restart) with contextual menu
- Kill method can be adjusted




14 / 18

Monitoring Gui + Browser




#### Monitoring Gui + Browser



- Browser to build histograms selections to display
- Search Function to refine selection.


Monitoring Gui + Browser



• List of histograms added from Browser




Monitoring Gui + Browser



Multiple canvas area available




Monitoring Gui + Browser



Multiple canvas area available



Monitoring Gui + Browser




- Multiple canvas area available
- Real ROOT histograms (Can be fitted, zoomed, etc.)



16 / 18

Monitoring Gui + Browser




- Multiple canvas area available
- Real ROOT histograms (Can be fitted, zoomed, etc.)
- Histograms descriptions and Quality



16 / 18


Monitoring Gui + Browser



- Multiple canvas area available
- Real ROOT histograms (Can be fitted, zoomed, etc.)
- Histograms descriptions and Quality



Monitoring Gui + Browser



- Multiple canvas area available
- Real ROOT histograms (Can be fitted, zoomed, etc.)
- Auto Update



# Conclusion and plans

Conclusions and plans

### Conclusion and plans

- Independent processes decoupled and linked using networking.
- Plugins (modules, data streaming) to configure and run the system.
- Tools for data feed in the system from the DAQ (event client interface)
- GUIs to control/monitor the system.
- Tests are OK but need numbers!
- II CSOFT release?

### Current work and plans

- Full implementation of SDHCAL DQM
- Combined ECAL test-beam -> combined ECAL-HCAL DQM as proof of concept
- EUDAQ binding (T. Coates, Sussex, UK)



# Backups

#### xdrstream and xdrlcio

xdrstream github: https://github.com/DQM4HEP/xdrstream
 Serialization with XDR format to read/write raw data into different devices (file, buffer, socket, user defined)



### Globally -> SIO implementation :

- without static API (SIO blockManager, SIO streamManager, etc)
- with different device implementation
   external standalone package

For the moment, only buffer implementation (xdrstream::BufferDevice)

- xdrlcio github: https://github.com/DQM4HEP/xdrlcio LCIO edm serialization using xdrstream.
  - •XdrLcio::writeEvent(const EVENT::LCEvent \*pLCEvent, xdrstream::IODevice \*const pDevice)

#### and

- •XdrLcio::readNextEvent(xdrstream::IODevice \*const pDevice)
- + many useful functions

